


Transepitelial Multi-Use®



Pilares Multi-Use®

Pilar Recto Multi-Use®

TITANIO 5 NP 3.5

30 Ncm



h1.5 mm
VCON10N1500



h2.5 mm
VCON10N2500



h3.5 mm
VCON10N3500



h4.5 mm
VCON10N4500



h5.5 mm
VCON10N5500



h6.5 mm
VCON10N6500

Pilar Angulado Multi-Use®

TITANIO 5 NP 3.5

30 Ncm



17° h2.5 mm
VCON10N2517



17° h3.5 mm
VCON10N3517



30° h3.5 mm
VCON10N3530



30° h4.5 mm
VCON10N4530

Cicatrización

Tapón de Cicatrización Multi-Use® (Tornillo integrado)



TITANIO 5 MU0102

Impresión

Coping de Impresión Multi-Use®

TITANIO 5



Cubeta Abierta
MU0211



Cubeta Cerrada
MU0202


Réplica Multi-Use®

ACERO INOX



MU03

Transepitelial Multi-Use®



Prótesis Atornillada

Calcinable
Multi-Use®



POM

MU0402


Pilar de Sobrecorado
Multi-Use® Recto



POM CR
CO

MU0502

Pilar Provisional
Multi-Use®



TITANIO
5

MU0602

PEEK

MU0602P

Tornillo Protésico Multi-Use®



MAX 15 Ncm



TITANIO
DLC

MU0905T



TITANIO
5

MU0905

Prótesis Atornillada Angulada Multi-Use®

Pilar de Sobrecorado Multi-Use® Angulado



1



2

17°

1 VCHA170001
2 MU0504



1



2

30°

1 VCHA300001
2 MU0504

Tornillo Protésico Multi-Use® Angulado



MAX 15 Ncm



TITANIO
5

MU0900-TLB



Corta
VSDTLB-1
Mediana
VSDTLB-2
Larga
VSDTLB-3

Llave recta Multi-Use®



Contra-Ángulo
VDMU-1



Carraca
VDMU-2

Multi-Use® Pilar VulkanLoc



MU11R00